



Molecular sieve vacuum swing adsorption for gaseous dark matter and rare-event detectors

IOP Joint APP, HEPP and NP 2024

R.R. Marcelo Gregorio, N.J.C. Spooner, F. Dastgiri, A. C. Ezeribe, G. Lane,
A.G. McLean, K. Miuchi, H. Ogawa

Queen Mary University of London
University of Sheffield

April 2024
Liverpool, UK



Introduction

MS VSA for
Gas-Based
DM detectors

Introduction

Directional
Detectors

Radon


Common
Pollutants

Molecular
Sieves

Gas System

Further Work

Summary



- Next-generation direct detection dark matter experiments increase in sensitivity, they approach an irreducible neutrino background
- If no dark matter signals are detected, this background establishes an ultimate discovery limit for these experiments
- Directional detectors offer a way to mitigate this issue



Directional Detectors

Extra Information

MS VSA for
Gas-Based
DM detectors

Introduction

Directional
Detectors

Extra Information

Gas TPCs
CYGNUS

Radon



Common
Pollutants

Molecular
Sieves

Gas System

Further Work

Summary





Directional Detectors

Gas TPCs

MS VSA for
Gas-Based
DM detectors

Introduction

Directional
Detectors

Extra Information

Gas TPCs

CYGNUS

Radon


Common
Pollutants

Molecular
Sieves

Gas System

Further Work

Summary



- Gas-based Time Projection Chambers (TPCs) are the most commonly employed for directional searches
- These detectors use low-pressure gas targets to produce long nuclear recoil tracks
- SF_6 has gained popularity in gas-based directional detectors due to its novel properties ^a

^aN.S. Phan et al 2017 JINST 12 P02012



Directional Detectors

CYGNUS Collaboration

MS VSA for
Gas-Based
DM detectors

Introduction

Directional
Detectors

Extra Information

Gas TPCs

CYGNUS

Radon


Common
Pollutants

Molecular
Sieves

Gas System


Further Work

Summary



CYGNUS-HD10
Strip micromegas
He:CF₄:X
40 L + 1 m³ R&D
detectors under
construction

CYGNUS-10
10 m³, GEMs + wires
He:SF₆
Boulby, UK
R&D ongoing on 1 m³



CYGNUS/INITIUM
GEMs + sCMOS + PMT
He:CF₄ (:SF₆)
LNGS, Italy
1 m³ demonstrator
funded towards 30 m³
experiment



CYGNUS-KM
1 m³, GEMs + 2D strips
SF₆/CF₄
Kamioka, Japan
R&D ongoing on 1 m³



CYGNUS-OZ
Stawell, Australia
GEMs + CCDs for gas studies
Small prototype under
development





Radon Contamination

MS VSA for Gas-Based DM detectors

Introduction

Directional Detectors

Radon

Contamination

Backgrounds


Common Pollutants

Molecular Sieves

Gas System

Further Work

Summary



- Radon is a naturally occurring radioactive noble gas
- Contamination in detectors originate intrinsically from detector material or plate out from atmospheric radon
- Radon contamination can act as a source of unwanted background noise from its decay and progeny decay



Radon

Backgrounds

MS VSA for
Gas-Based
DM detectors

Introduction

Directional
Detectors

Radon

Contamination

Backgrounds


Common
Pollutants

Molecular
Sieves


Gas System

Further Work


Summary



(a) Nuclear recoil



(b) Neutron emission



(c) Electron recoil



Common Pollutants

MS VSA for
Gas-Based
DM detectors

Introduction

Directional
Detectors

Radon

Common
Pollutants


Molecular
Sieves

Gas System

Further Work

Summary

GEM - O₂ Test - Detector Gain



- Common pollutants (O₂, N₂ and H₂O) can reduce the detector's amplification capabilities (Gain) by capturing electrons produced during interactions^a
- SF₆ is the most potent greenhouse gas, making the method of continuous flow and disposal problematic
- In industry gas filters are readily available to remove common pollutants

^aR Guida et al 2020 J. Phys.: Conf. Ser. 1498 012036



Molecular Sieves

Radon Removal

MS VSA for
Gas-Based
DM detectors

Introduction

Directional
Detectors

Radon

Common
Pollutants

Molecular
Sieves


Radon Removal

Low Background MS

Gas System

Further Work

Summary



- The Molecular Sieve was able to reduce the initial radon of $\sim 4 \text{ kBq/m}^3$ in SF₆ gas by 87%¹

¹R.R. Marcelo Gregorio et al., *Demonstration of radon removal from SF₆ using molecular sieves*, JINST 12 P09025 (2017)



Commercial Molecular Sieves are DIRTY!

- We want **High** Radon Filtration but **Low** Radon Emanation
- Commercial MS 5.4×10^{-3} Radon emanated per captured

Table 1. Radon filtration, intrinsic MS emanation and comparison parameter results for the NU-developed MS in granule and powdered form and NU-developed MS V2.

NU-developed MS	^{222}Rn Captured per kg (Bq kg^{-1})	^{222}Rn Emanated per kg (mBq kg^{-1})	^{222}Rn Emanated per ^{222}Rn Captured ($\times 10^{-3}$)
V1 (Granules)	35 ± 2	99 ± 23	2.8 ± 0.7
V1 (Powder)	330 ± 3	680 ± 30	2.1 ± 0.1
V2 (Powder)	254 ± 3	< 14.4	$< 5.7 \times 10^{-2}$

- Nihon University MS V2 emanates **99%** less per captured compared to commercial ²

²R.R. Marcelo Gregorio et al., *Test of low radioactive molecular sieves for radon filtration in SF6 gas-based rare-event physics experiments*, JINST 16 P06024 (2021)



Gas System Design

MS VSA for
Gas-Based
DM detectors

Introduction

Directional
Detectors

Radon

Common
Pollutants

Molecular
Sieves

Gas System

Design

Prototype Demo


Radon Reduction

Gain Conservation

Further Work

Summary

The dual MS column configuration ensures an MS filter is always available by allowing simultaneous filtration and regeneration





Gas System

Prototype Demonstration

MS VSA for
Gas-Based
DM detectors

Introduction

Directional
Detectors

Radon

Common
Pollutants

Molecular
Sieves

Gas System
Design

Prototype Demo

Radon Reduction


Gain Conservation

Further Work

Summary

To assess the performance of MS, a gas system prototype was constructed and applied to a lab-based TPC detector with Thick Gas Electron Multiplier (ThGEM)

- Radon reduction test
- Gas gain conservation test





Gas System

Radon Reduction

MS VSA for
Gas-Based
DM detectors

Introduction

Directional
Detectors

Radon

Common
Pollutants

Molecular
Sieves

Gas System

Design

Prototype Demo


Radon Reduction

Gain Conservation

Further Work

Summary

- The gas system prototype resulted in suppressing the radon activity to 0.8 ± 6.4 mBq
- The large error margin can be attributed to the background limits of radon measurement apparatus
- The intrinsic radon activity in the ThGEM-based TPC detector setup to **less than 7.2 mBq** at a 95% C.L





Gas System

Gain Conservation

MS VSA for
Gas-Based
DM detectors

Introduction

Directional
Detectors

Radon

Common
Pollutants

Molecular
Sieves

Gas System

Design



Prototype Demo

Radon Reduction

Gain Conservation

Further Work

Summary





Gas System

Gain Conservation

MS VSA for
Gas-Based
DM detectors

Introduction

Directional
Detectors

Radon

Common
Pollutants

Molecular
Sieves

Gas System

Design

Prototype Demo


Radon Reduction

Gain Conservation

Further Work

Summary

- Gain stayed above background until detector operation was terminated
- Appears to approach a steady state gain
- ~ 3 times less gas compared to conventional fresh gas flush systems





Further Work

MS VSA for
Gas-Based
DM detectors

Introduction

Directional
Detectors

Radon

Common
Pollutants

Molecular
Sieves

Gas System

Further Work

Summary

- Gas System and MS are currently being implemented to directional dark matter experiments (NEWAGE, Japan, CYGNO, Italy and CYGNUS-OZ, Australia)
- Awarded a Grant to build a Radon Emanation Facility in Australia (DREAMR)
- Next generation filters Metal Organic Framework (MOFs) research
- Investigation to filter application to other target material such as Xenon and Argon



Summary

MS VSA for
Gas-Based
DM detectors

Introduction

Directional
Detectors

Radon

Common
Pollutants


Molecular
Sieves

Gas System

Further Work

Summary

- Two independent problems arise in gas-based DM detectors due to contaminant gases
- Radon contamination can mimic genuine signals
- Common pollutants can reduce gain capabilities
- These are two separate issues can be addressed with appropriate low radioactive MS filters
- A prototype gas system was constructed and applied to a ThGEM-based TPC detector, demonstrating reduction of intrinsic radon contamination and conservation of detector gain





Thank You! Any Questions?

MS VSA for
Gas-Based
DM detectors

Introduction

Directional
Detectors

Radon

Common
Pollutants

Molecular
Sieves

Gas System

Further Work

Summary

Dr. Robert Renz Marcelo Gregorio, CPhys

✉ robert.gregorio@cern.ch

in [linkedin.com/in/rrmg](https://www.linkedin.com/in/rrmg)

@ [@renzmarcelo](https://www.instagram.com/renzmarcelo)

List of Relevant Papers in This Work:

- R.R. Marcelo Gregorio et al., JINST 12 P09025
- R.R. Marcelo Gregorio et al., JINST 16 P06024
- R.R. Marcelo Gregorio et al., JINST 19 P03012

